

Accu-Chek® FastClix

The world's first 1-Click lancing device with a drum

£2.45
for
102
Lancets*



Key points:

- No clinical waste (subject to following local guidance)
- Prime and release in 1-Click action - easy to teach
- 11 individual depth settings
- Clixmotion technology provides gentle lancing and helps reduce the risk of skin and nerve-end damage.
- Virtually pain free testing
- 6 preloaded lancets in a drum, no single lancets to see or touch
- Penetration range: 0.7 to 2.2 mm
- Prime & release of device: press once to lance
- Size of device (length x width): 119 x 19 mm
- AST capability: Yes (with AST cap, available on request)
- Needles are enclosed in a drum through loading and unloading to avoid adverse needle stick injuries.

Available with all Accu-Chek meters and insulin pumps, also available as a standalone device.



A starter pack contains everything a patient needs to get started straight away, making it easy to use.

Roche Diabetes Care products offer quality and accuracy from the market leading brand, with over 40 years' experience in diabetes care. We partner with healthcare customers using best in class implementation to create prescribing efficiencies that add value to diabetes care within health economies.



For more information
RocheDiabetes.co.uk
RocheDiabetes.ie

© 2020 Roche Diabetes Care Limited. All rights reserved.
ACCU-CHEK, ACCU-CHEK AVIVA EXPERT, ACCU-CHEK PERFORMA, ACCU-CHEK MOBILE and FASTCLIX are trademarks of Roche. All other trademarks or brand names are the property of their respective owners.

Roche Diabetes Care Limited, Charles Avenue, Burgess Hill, West Sussex, RH15 9RY. Company Registration Number: 9055599

Date of preparation: September 2020

For healthcare professional use only in the UK and Ireland.
Not for distribution.

www.accu-chek.co.uk www.accu-chek.ie

If you are interested in learning more about how Roche Diabetes Care can support with formulary implementation contact our dedicated support team on UK: **0800 0407221**

*Based on list price of £5.90 for 204 lancets on the UK Drug Tariff as of June 2020.

Safety advice: The Accu-Chek FastClix finger pricker is intended for patient self-monitoring by an individual person only. It must not be used to collect blood in a multipatient setting as it does not incorporate any features to guard against cross infection.

ACCU-CHEK®